


SELECTING CULTURE MEDIA

Biplate

Results: Gram-positive, Gram-negative bacteria, No growth or contaminated

Factor Agar
Growth of Gram-Positive
bacteria




MacConkey Agar
Growth of Gram-Negative
bacteria

Triplate

Results: Staphylococcus spp., Streptococcus spp., Gram-negative bacteria,
No growth or contaminated

Factor Agar
Growth of Gram-Positive
bacteria




MacConkey Agar
Growth of Gram-Negative
bacteria

Modified TKT (MTKT) Agar
Growth of *Streptococcus*

Quad plate

Results: Staphylococcus spp., Streptococcus spp., Gram-negative bacteria,
No growth, contaminated or others.

Blood Agar
Growth of most bacteria



MacConkey Agar
Growth of Gram-Negative
bacteria

Factor Agar
Growth of Gram-Positive
bacteria

Modified TKT (MTKT) Agar
Growth of *Streptococcus*